

# INSERT WEAR GUIDE

How the best perform



After Five® Right

After Five® Straight

After Five® Left

#3 Beaver tail

Swivel XT™

#1000 Triple Bend

#100 Thin

#10 Universal

Original Propriety Designs

Furcation Right

Furcation Straight


Furcation Left



Hu-Friedy's NEW Swivel XT -  
the thin insert that's tough  
on heavy calculus

### INSERT WEAR GUIDE INSTRUCTIONS

1. Locate the diagram of the insert to be measured by color.
2. Position point **A** to point **B** as shown.
3. Rotate insert until tip is flat against the card.
4. Evaluate tip condition, comparing length to red/blue bars.



Tip extends beyond blue line  
**Optimum Efficiency**  
(Scale Away)



Tip touches blue line  
**25% Efficiency Loss**  
(Reorder)



Tip touches red line  
**50% Efficiency Loss**  
(Discard)

To reorder inserts call  
1-800-Hu-Friedy or your  
authorized Hu-Friedy Dealer



©2011 Hu-Friedy Mfg. Co., LLC. All Rights Reserved. Hu-Friedy, the Hu-Friedy oval logo, After Five and Swivel XT are the trademarks of Hu-Friedy Mfg. Co., LLC. HF-4840-01H